

Welcome
to the world of
smart parenting


Why does
your baby need
a **BEMPU**
Bracelet ?

Why should I use
the BEMPU Bracelet
for my baby?

??



Second, cold stress often
goes unnoticed!



bempu



Instructions for use


BEMPU is a smart baby bracelet.

It alerts you to take action if your
baby is at risk of becoming
hypothermic.



Well first, now that your baby
is outside of your womb she
can have trouble maintaining
her body temperature by
herself and this leads to
'cold stress'

..& a precise and warm
body temperature is
essential for healthy
development.



Luckily, the BEMPU
bracelet alerts you
if your baby needs
warmth!


.. & a warm baby
grows to be a
healthy baby!!




How to use your BEMPU Bracelet




Step 1 - Pull out the tag to turn the device ON



The bracelet will play a welcome tune and blink light, indicating the device is **ready to use**.




Step 2 - Put the bracelet on baby's wrist.



Blue Light ● ⇒ Your baby's temperature is OK ✓


The bracelet shows a BLUE colour if your baby's temperature is normal.



Alarm goes ON
&
Light turns Orange ●

⇒ Your baby needs warmth !!


The bracelet **SOUNDS ALARM** and indicates an **ORANGE** colour if your baby's temperature is low and she needs warmth!



When light is Orange ●

Hold the baby directly against your skin to warm the baby. This is called **Kangaroo Care**.

Provide skin-to-skin contact !!




When light is Orange ●

Swaddle in a warm cloth

Also, swaddle your baby in a warm cloth and give her a hat, socks and mittens.

The bracelet will return to show a blue light when the baby becomes warm.




Light turns **Blue** ●

Regular kangaroo care is a proven way to promote thermal control, prevent infections, and encourage healthy growth.



for
24/7 1 month

Have your baby wear the BEMPU bracelet for 1 month and keep her protected from hypothermia.



Do's & Don't's


Regularly clean the BEMPU bracelet

1. Clean the bracelet by wiping it with a soft wet tissue.
2. All wipes which are suitable for baby use can be used for cleaning.
3. DO NOT sterilize the bracelet or put it in water.


Put the BEMPU bracelet snugly on the wrist



Snug fit but not loose



Not too tight



Disposing the bracelet

Dispose the device after using for 1 month

1. Stop using the device after 1 month.
2. After 1 month the device will constantly show a white light for 1 full day and then will go dark, indicating its end of use.
3. Dispose the device such that it is not in the reach of children.
4. DO NOT burn or crush the device for disposing.